

SCR Characteristics Apparatus with Regulated Power Supply

Objective:

Product Code:GE221

To plot Anode Current Vs Anode Cathode Voltage Characteristics. Gate Characteristics (Open Gate), Gate Characteristics (Gate connected).



Specifications:

- Instrument comprises of Two DC Regulated Power Supplies 0-30VDC/1 50 mA & 0-1 VDC/150mA
- Three Round Meters for Voltage & Current Measurement
- SCR 2P4M mounted behind the panel
- Connections of Supplies, Meters & SCR brought out at 4mm Sockets.