

Photo Transistor Characteristics Apparatus

Product Code :GE 233

Objective: To Plot A Graph, Current Vs Intensity OF Light At Different Value OF Voltages.



Specifications:

- Instrument comprises of DC Regulated Power Supply 0-3 VDC/150mA,
- Two Round Meters For Voltage & Current Measurements,
- Photo Transistor mounted in a Metal box
- Connections of Supplies & Photo Transistor brought out at 4mm Sockets.
- Lamp Holder with 60W & 50 cm Wooden plank are given with this apparatus.