

## Characteristics of Opto Electronic Devices

**Product Code :GE 234**

**Objective:** To Study V-I Characteristics of LDR, Photo Transistor, Photo Diode & LED.



### Specifications:

- Instrument Comprises of DC Regulated Power Supply 0-3 VDC/150mA,
- Two Round Meters for Voltage & Current Measurement,
- LDR, Photo Transistor, Photo Diode & LED mounted in a Metal box,
- Connections of Supplies, Meters & Components brought out at 4mm Sockets.
- Lamp Holder with 60 W Bulb & 50cm Wooden Plank are given with this apparatus.