



GURU KRIPA ENTERPRISES

AN ISO 9001:2015 & NSIC CERTIFIED COMPANY

Sur No - 668, Plot No - 268 , Fatevadi, Sarkhej , Ahmedabad,Gujarat - 382210

Email: info@gkelabs.com, www.gkelabs.com, Mobile::9106316361



GEM PORTAL
REGISTERED

Conversion of Galvanometer into a Voltmeter/Ammeter

Product Code: GE247

Objective To Calculate the value of Series Resistance for Voltmeter & Parallel Shunt for Current Meter & to Convert Galvanometer to Voltmeter & Ammeter.



Features : Instrument comprises of DC Regulated power supply of 0-5V. Three round moving coil Meters (One Voltmeter, One Milliammeter & One Galvanometer). Three shunts, Three Resistances are mounted on front panel.

