



GURU KRIPA ENTERPRISES

AN ISO 9001:2015 & NSIC CERTIFIED COMPANY

Sur No - 668, Plot No - 268 , Fatevadi, Sarkhej , Ahmedabad,Gujarat - 382210

Email: info@gkelabs.com, www.gkelabs.com, Mobile::9106316361



GEM PORTAL
REGISTERED

RAM Circuit using IC 7489

Product Code:GE392

Objective: To Study and verify the truth table of RAM (Random Access Memory) circuit.



Features: Instrument comprises of DC Regulated Power Supply 5VDC/1 50mA, 10 logic inputs selectable using SPOT switches, four LED output indicators, IC 7489 placed on the front panel & connections brought out at sockets.

