

TRANSMISSION LINE TRAINER

Productcode: GKE446



Features:

- Transmission Line: Coaxial cable 100m (25meterx4)
- Impedance Matching Resistances; 0100 Ohm, 2 Nos. •
- Test Generators: Sine wave 40KHz-400KHz (Low Range)
- Square wave 400KHz-4MHz

LIST OF EXPERIMENTS

- Measurement with Matched, Short & Open end of the line.
- Measurement of line properties.
- Measurement of line attenuation
- Frequency characteristics of line
- Input impedance of the line
- Measurement of stationary waves
- Phase shift along the line
- (Fault localization within the line
- Line in pulsed condition.