

Function Generator 0.1 Hz to 1 MHz Frequency counter

Product Code: GE514



- Waveshapes: Sine, Square & Triangular Wave outputs at 4mm sockets/terminals.
- Frequency Range: 0.1 Hz-1 MHz adjustable using fine & coarse controls.
- Accuracy: $\pm 3\%$ on all the ranges
- Amplitude for Sine Wave output: 30V peak to peak
- Amplitude for Square Wave output: 15V peak to peak
- Output impedance for Sine Wave: App. 600 Ohms
- Output impedance for Square Wave: App. 400 Ohms
- Display: Frequency