

Velocity of Ultrasonic waves in Kerosene Liquid

Product Code: GE 464



Ultrasonic Interferometer, a simple and direct device which are designed and fabricated for determining the ultrasonic velocity in liquids with a high degree of accuracy in the field of laboratory purpose in engineering colleges. Besides. the company is one of the most reputed Ultrasonic Interferometer suppliers in India and abroad.

Ultrasonic Interferometer of Liquid Consists of:-

- RF Oscillator 3-8 Mhz
- Specially Design Long Tube Spectrometer
- Crystal Tank
- Sodium Lamp
- Sodium Lamp Transformer
- Wooden Box for sodium lamp.