

Logic Gates using Universal gates (NAND & NOR Gates)

Product code :GE 377

Objective : To verify the truth table of AND gate, OR gate, NOR gate, NAND gate & NOT gate using NAND & NOR gate..



Features: Instrument comprises of DC Regulated Power Supply 5V/150mA, 2 logic 1 & 2 logic 0 inputs are provided on sockets, 2 Red LED output indicators, circuit diagram printed for 4 NAND & 4 NOR gates & their respective ICs placed inside the cabinet & connections brought out at sockets.