


Study Of 4 Bit Adder And Subtractor

Product Code: GE 378

Objective: To study

- Half Adder Using TTL Gates, Full Adder Using TTL Gates, Half Subtractor Using TTL Gates, Full Subtractor Using TTL Gates, 4 Bit Binary Adder Using 7483 IC , 4 Bit Binary Parallel Adder Using 7483 IC..



Features: Instrument comprises of DC regulated power supply 5V/150mA, 9 logic inputs selectable using SPDT switches, 5 output led indicators, circuit diagram printed & connections brought out at socket on the front panel.