

## IV-E237 To Determine The Mechanical Equivalent Of Heat By Joule's Calorimeter



### Specifications :

- It consists of a copper vessel of 3"x2" with outer vessel as a superior quality teak wood box of size 9.5 x 9.5 x 11.5 cm fitted with bakelite lid with hole for thermometer and stirrer.

### Scope of Supply :

- Joule's Calorimeter Trainer Kit
- Thermometer
- Joule's Calorimeter
- Rheostat
- Joule's Calorimeter Trainer Kit

### Specifications :

- AC Power Supply 0-12V, 3A
- AC Voltmeter 20V
- AC Ammeter 3A
- Provision to Attach Rheostat
- Output Termination to Calorimeter

### Optional Requirement :

- Digital Balance