

IV-E383 Digital Gauss Meter



Gauss Meter is an instrument to find the Magnetic Field of Any Object

Experiment :

- To find the Magnetic Flux (Electromagnetic Field of any Object)

Digital Gauss Meter :

- Range: 0-2 K Gauss & 0-20 K Gauss,
- Resolution: 1 Gauss at 0-2 K Gauss Range,
- Accuracy: $\pm 0.5\%$,
- Display: 3 ½ Digit, 7 Segment LED,
- Power: 220V, 50Hz
- Special Features: Indicates the direction of the Magnetic field.
- Gauss Probe

Salient Features :

- Front panel built with high class insulated Printed Circuit Board sheet with well printed circuits and symbols.
- Instruction manual.
- Connections are brought out through 3Pin Connector
- The trainer is housed in ABS Plastic and Metal Cabinet.
- Size of the trainer set 12"x4.5"