

## IV-E653 Travelling Microscope



### About Spectrometer:

- The travelling microscope consists of a cast iron base with machined vee-top surface and is fitted with 3 levelling screws, one being fixed.
- The microscope can travel horizontally 18cm and vertically 15cm with the help of slides.
- The slow motion knobs are provided for taking accurate readings.
- Vernier readings 0.01 mm.
- The Microscope tube consists of an eyepiece 10X with 50mm or 75 mm objective, for holding objects horizontal stage, made of milky cenolite sheet provided on the base.
- Packed in beautiful wooden box.