

IV-E1040 Photo Detector With Sensor



Specifications:

- Operated on 230V, 50Hz
 - Inbuilt Power Supply +5V, 150mA
 - Digital Readout in Amperes at different Ranges
 - Connections are brought out for Photo Detector
 - LDR / Photo Transistor / Photo Voltaic Cell is Used as Photo Sensor
 - Multi Range Meter
 - Rotary Switch for Selection of Ranges
- 1.200uA
 - 2.2mA
 - 3.200mA
 - 4.2000mA