

IV-E846 To Measure The Diameter Of A Wire With The Help Of Micrometer Screw Gauge

Specification:

- This economical instrument is sufficiently accurate for beginners and for most elementary laboratory measurements.
- Superior quality
- Semi circular Welch pattern
- Dull nickel finish
- Oxidized threaded rod
- Individually packed in card board box
- With detailed working manual.



Size:

L.C.	0.01mm	0.001inch
A.	0-25mm	0-1 inch
B	25-50mm	1-2 inch
C	50-75mm	2-3 inch
D	75-100mm	3-4 inch
E	100mm	4-5 inch

Apparatus Required:

- Micrometer Screw Guage of any type
- Small piece of wire
- Coin & lead shot.