

## IV-E845 To Measure The Diameter Of A Small Spherical Body With The Help Of Vernier Caliper



### Specification:

- Quality steel chrome plated.
- The metric scale is graduated in millimeter for 12 cm
- The 10-part Vernier reads is 0.1 mm.
- The English scale is graduated to 1/16 inch for 5 inches its 8-part Vernier reads to 1/28 inch with adjusting wheel and depth gauge.
- Elegantly packed in card board box, with detailed working manual.
- Least Count 0.01 cm.
- Measuring range 125mm/5inch, Graduation 0.1mm/1/28inch

### Size:

8"	0.01cm
10"	0.01cm
12"	0.01cm

### Apparatus Required:

- Vernier Calliper
- Glass Beaker
- Cylindrical Body
- Pendulum Bob.