

## IV-75 DIGITAL LEAK TEST APPARATUS

### PRODUCT FEATURES

- Programmable Time With Digital Display
- Auto Stop After Completion of Set Time With Alarm
- Precise Control of Vacuum
- Heavy Duty Oil-free Vacuum Pump
- 1mm Thick SS-304 cabinet

### ACCESSORIES

- Instruction manual
- 12 inch vacuum desiccator
- Connecting pipes



### Technical specification

<b>Type</b>	GMP Model, In-built Test Unit
<b>Time Control &amp; Range</b>	Digital, Programmable, 0.01 sec to 99 min 59 sec
<b>Test Vacuum Range</b>	Adjustable up-to 650 mm Hg
<b>Vacuum Display</b>	2.5" SS Vacuum Gauge
<b>Vacuum Control &amp; Release</b>	By Single Front Panel Knob
<b>Audio Alarm at completion</b>	Yes
<b>Vacuum Pump</b>	Oil free, MV 27" +-5% , 1/4HP Motor
<b>Test Unit Size</b>	200mm, 250mm & 300mm (any one) (Specify in Order)
<b>Power</b>	230V±10%AC-50hz
<b>Dimension</b>	600x450x200mm (approx)
<b>Weight</b>	8kg (approx)