

IV-770 TOUCH SCREEN DIGITAL MOISTURE BALANCE

PRODUCT FEATURES

- Long life IR based heating
- Highly stable and accurate
- Colored touch screen
- Fast heating

ACCESSORIES

- Touch Pen
- Sample Tray
- Main lead
- Instruction manual
- Weight box



IR SOURCE

IR SOURCE



Technical specification

Max. Capacity	120 g
Accuracy	0.005g
Readability	0.01%
Weighing sensor	HBM
Temperature range	40-230 C
Heating source	High efficiency IR source
Heating way	Gentle, standard, fast
Halt mode	Automatic, manual, timed
Result display	Test temperature, test time, moisture content, solid content, dry weight, dynamic curve, moisture regain, solid content, weight loss, dry weight, test method, sample name
Tray material \ Tray diameter	Stainless \ Φ90mm
Interface	Standard RS232 communication interlace, which can be connected to printers, PCs and other peripheral devices, and USB can be customized
User management	Password protection, multiple users, auto-lock
Display Screen	5"HD touch screen
Product size	370x215x195mm
Weight	5kg approx.