

IV-760 TOUCH SCREEN DIGITAL MOISTURE BALANCE

PRODUCT FEATURES

- Halogen based heating
- Highly stable and accurate
- Colored touch screen
- Fast heating

ACCESSORIES

- Touch Pen
- Sample Tray
- Main lead
- Instruction manual
- Weight box

Technical specification

HALOGEN SOURCE



Max. Capacity	120g
Accuracy	0.005g
Readability	0.01%
Weighing sensor	HBM
Temperature range	40-230
Heating source \ way	High efficiency IR source \ Gentle, standard, fast
Halt mode	Automatic, manual, timed
Result display	Test temperature, test time, moisture content, solid content, dry weight, dynamic curve, moisture regain, solid content, weight loss, dry weight, test method, sample name
Tray material \ diameter	Stainless \ Φ 90mm
Interface	Standard RS232 communication interlace, which can be connected to printers, PCs and other peripheral devices, and USB can be customized
User management	Password protection, multiple users, auto-lock
Display Screen	5"HD touch screen
Product size	370x215x195mm
Package size	445x300x330mm
Weight	5kg approx.