

## IV-730 DIGITAL AUTOMATIC VISCOMETER

### PRODUCT FEATURES

- Stepping Motor for Accurate, reliable operation
- Direct readout of all measurement parameters
- Auto Ranging
- Time Function for measurement
- Accurate calibration

### ACCESSORIES

- Spindle 4 Nos
- Instrument stand 1Nos
- RS-232 cable 1Nos
- Instruction manual
- Power adapter



### Technical specification

<b>Measuring range</b>	10-2,000,000 C.P
<b>Display Resolution</b>	0.01 C.P
<b>Measurement accuracy</b>	±5% Newtonian Fluid)
<b>Speed</b>	0.3 0.6 1.5 3 6 12 30 60 RPM/min
<b>Outline size</b>	95\*130\*155 Excluding base Net weight 2kg Excluding base
<b>No of standard spandel</b>	#1, #2, #3, #4 is standard (#0 is optional)
<b>Display</b>	LCD display with backlit
<b>Power Supply</b>	230VAC ±10%, 50Hz
<b>Dimension</b>	400x350x360mm approx
<b>Weight</b>	9kg approx