

IV-151 DIGITAL TDS METER

INBUILT STIRRER

PRODUCT FEATURES

- Inbuilt stirrer
- Bright LED display
- Stable and accurate result
- Easy to operate

ACCESSORIES

- TDS cell 1.0
- Rotational TDS cell stand
- Magnetic capsule
- Instruction manual
- Dust cover



Technical specification

Display	3½ Digit LED
Measurement	TDS
Cond. Ranges	Five Ranges
Range Selection	Automatic 5 Ranges
Temp. Range	±5 to 110°C
Cond.	±0.5% FS ± 1DIGIT
Temp	±0.1°C ± 1DIGIT
Stirrer	Inbuilt Stirrer (With Speed Variable)
Temperature Compensation	Auto: 0 to 50 °C Manual: 0 to 50 °C
Cell Constant	Adjustable on Digital Display
Measuring Cell	Platinum DIP Type
Cell Const	Auto
Accessories	TDS Cell, Magnetic Capsule, Operation Manual, Dust Cover, Cell Stand,