

## IV-26 DIGITAL CONDUCTIVITY METER

### PRODUCT FEATURES

- Bright LED display
- Automatic temperature compensation
- Stable & accurate result
- Compact design

### ACCESSORIES

- Conductivity cell 1.0
- Temperature probe
- Rotational Conductivity cell stand
- Instruction manual
- Dust cover



### Technical specification

<b>Display</b>	3½ Digit LED
<b>Measurement</b>	Conductivity
<b>Cond. Ranges</b>	Five Ranges
<b>Range Selection</b>	Automatic in 5 Ranges
<b>Accuracy</b>	±0.5% FS ± 1DIGIT
<b>Temperature Compensation</b>	Auto 0 to 50°C
<b>Cell Constant</b>	Adjustable on Digital Display
<b>Measuring Cell</b>	Platinum DIP Type
<b>Cell Const. Measurement</b>	Auto
<b>Power</b>	230V±10% AC
<b>Dimensions</b>	270x202x80 mm (LXBXH)
<b>Weight</b>	0.8 Kg approx.