

IV-16 DIGITAL CONDUCTIVITY METER

PRODUCT FEATURES

- Bright LED display
- Manual temperature compensation
- Stable & accurate result
- Compact design

ACCESSORIES

- Conductivity cell 1.0
- Temperature probe
- Rotational Conductivity cell stand
- Instruction manual
- Dust cover



Technical specification

Display	3½ Digit LED
Measurement	Conductivity
Cond. Ranges	Five Ranges
Range Selection	Automatic in 5 Ranges
Accuracy	±0.5% FS ± 1DIGIT
Temperature Compensation	Manual 0 to 50°C
Cell Constant	Adjustable on Digital Display
Measuring Cell	Platinum DIP Type
Cell Const. Measurement	Auto
Power	230V±10% AC
Dimensions	270x202x80 mm (LXBXH)
Weight	0.8 Kg approx.