

IV-22 DIGITAL POTENTIO METER

INBUILT STIRRER

PRODUCT FEATURES

- Inbuilt stirrer
- 3½ Digit LED display
- Highly accurate and stable
- Compact design

ACCESSORIES

- Calomel electrode
- Glass electrode
- Platinum electrode (optional)
- Silver electrode (optional)
- Magnetic capsule
- Rotational electrode stand



Technical specification

Display	3½ digit bright red LED display
Range	0 to ± 1999 mV
Resolution	1 mV
Accuracy	$\pm 1\text{mV} \pm 1\text{digit}$
Operating Temperature	10°C to 50°C
Inbuilt Stirrer	Available
Power Supply	230V $\pm 10\%$ AC, 50 Hz
Dimensions	270 X 248 X 101 mm (LxBxH)
Weight	2 kg