

IV-32 DIGITAL POTENTIO METER

PRODUCT FEATURES

- Highly stable & accurate
- 3½ digit LED display
- Rugged & economical
- Compact design

ACCESSORIES

- Calomel electrode
- Glass electrode
- Platinum electrode (optional)
- Silver electrode (optional)
- Rotational electrode stand



Technical specification

Range	0 - 1999mV
Resolution	1 mV
Accuracy	± 1 mV
Input Impedance	1012 Ohms.
Operating Temperature	10°C to 50°C
Display	3½ Digit 0.5 seven segment LED Indication with auto polarity indication
Power	230V ± 10 % AC, 50Hz
Dimensions	76x275x170 mm
Weight	2 kg. (approx)