

IV-34 DIGITAL TURBIDITY METERS

PRODUCT FEATURES

- Display-3½ digit LED
- 30mm flat bottom Test Tube
- Ranges 0-200 NTU
- Highly Stable & Accurate

ACCESSORIES

- Test Tube 30mm dia - 1 Set
- Instruction Manual - 1 No.
- Dust Cover - 1 No.



Technical specification

Display	3½ Digit LED
Detector	Photocell/Photodiode
Range	0 to 200 NTU in 2 ranges
Resolution	0.1 NTU
ACCURACY	± 2% FS ± 1digit
Sample System	60 x 25mm Clear Glass Test Tube
Light Source	LED
Power	230V ± 10%, 50Hz. A
Dimensions	280 x 175 x 110 mm
Weight	3kg approx