

IV-181 DISSOLVED OXYGEN METER

PRODUCT FEATURES

- Inbuilt stirrer
- Bright LED display
- Compact design
- Accurate and stable result

ACCESSORIES

- Gold/silver D.O Probe
- D.O membrane
- Rotational D.O probe stand
- Magnetic capsule
- Dust cover
- Instruction manual

Technical specification

INBUILT STIRRER



Range	0 to 20.0 ppm
Resolution	0.1 ppm
Accuracy	0.1ppm \pm 1 count
Temp. Compensation	0°C to 50°C manual
Compensation	\pm 3%
D.O. Sensor	amperometric gold/silver membrane type.
Stirrer	Inbuilt stirrer with speed control
Display	3½ digit LED display
Power	230V A.C 50 Hz
Weight	2 Kg