

IV-1261 MICROPLATE ELISA WASHER

PRODUCT FEATURES

- Large LCD display, intelligent software, easy to operate
- Strong washing function, bottom washing and two positions aspiration make completed washing effect
- 8-channel or 12 channel optional
- Needle positions adjustable and up to 10 different needle positions pre-settable
- Pressure, volume & time adjustable, washing times, soap time pre-settable Pressure, volume & time adjustable, washing times, soap time pre-settable, whole plate or single strip reading
- Self calibration and trouble alarm function
- Large memory to store up to 100 program

ACCESSORIES

- 4 wash bottle (2500ml)
- Connecting tube set
- Cleaning wire
- Main lead
- Instruction manual



Technical specification

| | |
|--------------------------|---|
| Plate | 48/96 well and trips (U.V. Flat bottom) |
| Residual Volume | ≤2µl |
| Dispensing volume | Adjustable from 50 to 900 µl |
| Wash cycles | 1-16 optional |
| Mode of wash | From 1 to 12, row and rank, selectable |
| Soak time | Up to 990 seconds adjustable |
| Shake time | Up to 990 seconds adjustable |
| Display | Large LCD |
| Channel | 3 channels for washing liquid, 1 channel for waste liquid |
| Memory | up to 100 programs can be stored |
| Interface | RS232 |
| Dimension | 450(L)mm x 310(W)mm x 170 (H)mm |
| Net Weight | 7.5kg |
| Power | AC220V±10%, 50/60Hz; AC110V±10%, 50/60Hz |