


IV-83 MOTORLESS MAGNETIC STIRRERS

PRODUCT FEATURES


- Digital display provided easy reading of the speed and time.
- The most advance magnetic coil technology for wear-free drive without moving parts.
- Automatic reverse rotation (Pulse function) every 30 second provided excellent mixing performance.
- Max stirring quantity (H₂O) of 800 ml.
- Speed range of 15 to 1500 rpm.
- Slim design with unique flat surface design.
- Water spill resistant design - IP 65.

ACCESSORIES

- Stirrer bar 2 piece, 9V 1.5 power adaptor, operation manual



Technical specification

Image	
Product Model	IV-83
Speed Limit	15-1500 RPM
Stirring Position	1
Increment and decrements	50 RPM
Display	Digital LED
Working Plate Material	SS Plate
Setting Dia.	105 mm
Maximum Stirring Capacity	1000 ml H ₂ O
Recommended Stirrer Bar Size	25x8mm PTFE Coated
Pulse Mode \ Timer	30 second auto reverse direction \ 99 Minuits and Infinite Mode
Last run Memory Save	Yes
Protection Class	Ip65
Permissible Ambient Temp \ relative Humidity	5 °C to 40 °C \ 80%
Voltage \ Frequency	100-240 VAC \ 50/60 Hz